# SciCan STATIM G4

Cassette Autoclaves



# **STERILIZATION AND DRYNESS**

at speeds faster than conventional chambered autoclaves



### **Fast Cycle Times**

STATIM 2000 G4 unwrapped cycle times are the world's fastest at 6 minutes. Wrapped cycles, including drying, can be completed in 35 minutes.

STATIM 5000 G4 can complete unwrapped cycles in as little as 9 minutes. Wrapped cycles, including drying, can be completed in 40 minutes.



### **G4** Technology

G4 equipment can connect and send information to your smart devices so you can have cycle data, instructions, and maintenance notifications at your fingertips.



#### Gentle

The fresh water feature guarantees fresh water with every cycle. State of the art steam technology is gentle on both solid and hollow instruments.



#### Reliable

For over 25 years, STATIM's unmatched speed, reliability, and dependability have made it one of the world's best selling cassette autoclaves.



# Capacity

STATIM G4 autoclaves not only offer fast sterilization times but also enough capacity to meet the needs of a busy practice.

STATIM 2000 G4 capacity allows for up to 5 pouches and the STATIM 5000 G4 allows for up to 12 pouches.



STATIM 2000 G4 Specifications p/n G4-121101							
<b>Unit size (L x W x H):</b> 19.5″ x 16.25″ x 6″ / 49.5 cm x 41.5 cm x 15 cm		Unwrapped	Wrapped	Rubber / Plastic	Air Drying		
Cassette internal dimensions: 11" x 7" x 1.4" / 28 cm x 18 cm x 3.5 cm Minimum shelf height for proper leveling: 8.5" / 22 cm Weight without water: 46 lbs / 21 kg	Sterilization Time	3.5 minutes	10 minutes	15 minutes	Air drying only		
	Sterilization Temp	135° C	135º C	121° C	Air drying only		
	Total Cycle Time*	6 minutes	14 minutes	20 minutes	Air drying only		

**Power consumption:** 115 V, 60 Hz, 11A

\*Cycle time is dependent on instrument load and does not include drying.

#### STATIM 5000 G4 Specifications p/n G4-201103

	Pr				
<b>Unit size (L x W x H):</b> 23.6" x 16.1" x 7.5" / 60 cm x 41 cm x 19 cm		Unwrapped	Wrapped	Rubber / Plastic	Heavy Duty Unwrapped
<b>Cassette internal dimensions:</b> 15" x 7" x 3" / 38 cm x 18 cm x 7.5 cm	Sterilization Time	3.5 minutes	6 minutes	35 minutes	6 minutes
Minimum shelf height for proper leveling: 10" / 26 cm	Sterilization Temp	132º C	132º C	121º C	132° C
<b>Weight without water:</b> 73 lbs / 33 kg	Total Cycle Time*	9 minutes	17.5 minutes	45 minutes	12 minutes

\*Cycle time is dependent on instrument load and does not include drying.

**Power consumption:** 115 V, 60 Hz, 11A



**VistaPure** Water Filtration

Provides two grades of high-quality water for autoclaves, dental bottles, ultrasonic cleaners and the rinse cycle of instrument washers.

#### V3000



**ENSURE** Monitoring for steam sterilization

An easy to use, in-office monitoring system for steam sterilization.

SCI-BISK-120US – Ensure BI Starter Kit (25 BI's, Incubator & Log Book)



## VistaPure AutoFill

Replenish any new or older-model STATIM and BRAVO (classic unit) autoclave with distilled-quality water automatically, enabling truly hands-free convenience.

S3020 – AutoFill for STATIM S3030 – AutoFill for BRAVO

**STATMAT** 

Protect your STATIM

01-1128675 - For STATIM 2000

01-1128685 - For STATIM 5000

Securely seat cassettes atop your STATIM

with a safe and hygienic silicone mat.



#### **VistaCool** Direct-to-Drain System

Protects cabinets and drains—and saves valuable staff time—by cooling autoclave wastewater and sending it directly down the drain automatically.

V7501 – For One Autoclave V7502 – For Two Autoclaves V7601 – For One Autoclave

#### **Coltene International Dental Group**

Manufactured for: Dent4You AG Bahnhofstrasse 2 CH-9435 Heerbrugg

Distributed by: SciCan Ltd. 1440 Don Mills Road Toronto, Ontario M3B 3P9 Canada Distributed by: Coltène/Whaledent Inc. 235 Ascot Parkway Cuyahoga Falls, Ohio 44223 US

VistaPure and VistaCool are manufactured by: Crosstex International, Inc. 1733 W. Parkside Lane Phoenix, AZ 85027 T 800 531-2221 VistaResearchGroup.com

#### www.scican.com

SD-531 US EN R2

STATIM is a registered trademark and ENSURE, STATMAT, and BRAVO are trademark of SciCan Ltd., used under license by Dent4You AG. All Vista product names are trademarks of HuFriedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

